

Irish Standard I.S. EN 17322:2020&LC:2020

Environmental Solid Matrices -Determination of polychlorinated biphenyls (PCB) by gas chromatography - mass selective detection (GC-MS) or electroncapture detection (GC-ECD)

© CEN 2020 No copying without NSAI permission except as permitted by copyright law.

I.S. EN 17322:2020&LC:2020

Incorporating amendments/corrigenda/National Annexes issued since pub	lication:
---	-----------

The National Standards Authority of Ireland (NSAI) produces the following categories of formal documents:

I.S. xxx: Irish Standard — national specification based on the consensus of an expert panel and subject to public consultation.

S.R.~xxx: Standard~Recommendation-recommendation~based~on~the~consensus~of~an~expert~panel~and~subject~to~public~consultation.

SWiFT~xxx: A~rapidly~developed~recommendatory~document~based~on~the~consensus~of~the~participants~of~an~NSAI~workshop.

This document replaces/revises/consolidates the NSAI adoption of the document(s) indicated on the CEN/CENELEC cover/Foreword and the following National document(s):

NOTE: The date of any NSAI previous adoption may not match the date of its original CEN/CENELEC document.

This document is based on: Published:

This document was published
under the authority of the NSAI
and comes into effect on:

13.030.01
13.030.10
2020-10-16
13.030.20
13.080.10

NOTE: If blank see CEN/CENELEC cover page

NSAI T +353 1 807 3800 Sales:

 1 Swift Square,
 F +353 1 807 3838
 T +353 1 857 6730

 Northwood, Santry
 E standards@nsai.ie
 F +353 1 857 6729

 Dublin 9
 W NSAI.ie
 W standards.ie

Údarás um Chaighdeáin Náisiúnta na hÉireann

This is a free page sample. Access the full version online.

National Foreword

I.S. EN 17322:2020&LC:2020 is the adopted Irish version of the European Document EN 17322:2020, Environmental Solid Matrices - Determination of polychlorinated biphenyls (PCB) by gas chromatography - mass selective detection (GC-MS) or electron-capture detection (GC-ECD)

This document does not purport to include all the necessary provisions of a contract. Users are responsible for its correct application.

For relationships with other publications refer to the NSAI web store.

Compliance with this document does not of itself confer immunity from legal obligations.

In line with international standards practice the decimal point is shown as a comma (,) throughout this document.

This is a free page sample. Access the full version online.

This page is intentionally left blank



Correction Notice

EN 17322:2020

Reference:

Title:	Environmental Solid Matrices - Determination of polychlorinated biphenyls (PCB) by gas chromatography - mass selective detection (GC-MS) or electron-capture detection (GC-ECD)					
Work Item:	00444167					
	Brussels, 2020-09-09					
Please include the following minor editorial correction(s) in the document related to:						
the following ☑ Engli ☐ Frend ☐ Gern	ch					
PQ/U						

It has been brought to our attention that this document, issued on 2020-07-29, requires modification.

On the title page, in the European foreword and in the Introduction, the reference to EN 16167:2018 has been changed to EN 16167:2018+AC:2019.

Please find enclosed the updated English version.

We apologise for any inconvenience this may cause.

This is a free page sample. Access the full version online. I.S. EN 17322:2020&LC:2020

This page is intentionally left BLANK.

EUROPEAN STANDARD

EN 17322

NORME EUROPÉENNE

EUROPÄISCHE NORM

July 2020

ICS 13.030.01; 13.030.10; 13.030.20; 13.080.10

Supersedes EN 15308:2016, EN 16167:2018+AC:2019

English Version

Environmental Solid Matrices - Determination of polychlorinated biphenyls (PCB) by gas chromatography - mass selective detection (GC-MS) or electron-capture detection (GC-ECD)

Matrices solides environnementales - Dosage des polychlorobiphényles (PCB) par chromatographie en phase gazeuse-spectrométrie de masse (CG-SM) ou chromatographie en phase gazeuse avec détection par capture d'électrons (CG-ECD) Feststoffe in der Umwelt - Bestimmung von polychlorierten Biphenylen (PCB) mittels Gaschromatographie und massenspektrometrischer Detektion (GC-MS) oder Elektronen-Einfang-Detektion (GC-ECD)

This European Standard was approved by CEN on 15 June 2020.

This European Standard was corrected and reissued by the CEN-CENELEC Management Centre on 9 September 2020.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN-CENELEC Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION COMITÉ EUROPÉEN DE NORMALISATION EUROPÄISCHES KOMITEE FÜR NORMUNG

CEN-CENELEC Management Centre: Rue de la Science 23, B-1040 Brussels

EN 17322:2020 (E)

Con	tents	Page
Euroj	pean foreword	3
Intro	duction	4
1	Scope	5
2	Normative references	5
3	Terms and definitions	6
4	Principle	
5	Interferences	
5.1	Interference with sampling and extraction	
5.2	Interference with GC	
6	Safety remarks	8
7	Reagents	9
7.1	General	
7.2	Reagents for extraction	
7.3	Reagents for clean-up	
7.4	Gas chromatographic analysis	
7.5	Standards	
7.6	Preparation of standard solutions	
8	Apparatus	
8.1	Extraction and clean-up procedures	
8.2	Gas chromatograph	16
9	Sample storage and preservation	16
9.1	Sample storage	16
9.2	Sample pre-treatment	16
10	Procedure	17
10.1	Blank test	
10.2	Extraction	
10.3	Concentration	20
10.4	Clean-up of the extract	20
10.5	Addition of the injection standard	
10.6	Gas chromatographic analysis (GC)	
10.7	Mass spectrometry (MS)	
10.8	Electron capture detection (ECD)	
11	Performance characteristics	30
12	Precision	30
13	Test report	30
Anne	x A (informative) Repeatability and reproducibility data	31
A.1	Materials used in the inter-laboratory comparison study	31
A.2	Inter-laboratory comparison results	32
Anne	x B (informative) Examples for gas chromatographic conditions and retention times	
	of PCBs	36
Anne	x C (informative) Calculation method for the estimation of total PCB content	37
Biblio	ography	45

EN 17322:2020 (E)

European foreword

This document (EN 17322:2020) has been prepared by Technical Committee CEN/TC 444 "Environmental characterization of solid matrices", the secretariat of which is held by NEN.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by January 2021, and conflicting national standards shall be withdrawn at the latest by January 2021.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN shall not be held responsible for identifying any or all such patent rights.

This document supersedes EN 16167:2018+AC:2019 and EN 15308:2016.

This document is the result of the merging of EN 16167:2018+AC:2019 and EN 15308:2016, with minor technical modifications.

According to the CEN-CENELEC Internal Regulations, the national standards organisations of the following countries are bound to implement this European Standard: Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Republic of North Macedonia, Romania, Serbia, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom.

EN 17322:2020 (E)

Introduction

Polychlorinated biphenyls (PCB) have been widely used as additives in industrial applications where chemical stability has been required. This stability on the other hand creates environmental problems when PCB are eventually released into the environment. Since some of these PCB compounds are highly toxic, their presence in the environment (air, water, soil, sediment and waste) is regularly monitored and controlled. At present determination of PCB is carried out in these matrices in most of the routine laboratories following the preceding steps for sampling, pre-treatment, extraction and clean-up, by measurement of specific PCB by means of gas chromatography in combination with mass spectrometric detection (GC-MS) or gas chromatography with electron capture detector (GC-ECD).

This document was developed by merging of EN 16167:2018+AC:2019, initially elaborated as a CEN Technical Specification in the European project 'HORIZONTAL' and validated by CEN/TC 400 with the support of BAM, with EN 15308, published by CEN/TC 292.

Considered the different matrices and possible interfering compounds, this document does not contain one single possible way of working. Several choices are possible, in particular relating to clean-up. Detection with both MS-detection and ECD-detection is possible. Two different extraction procedures are described and 9 clean-up procedures. The use of internal and injection standards is described in order to have an internal check on choice of the extraction and clean-up procedure. The method is as far as possible in agreement with the method described for PAH (EN 16181:2018 and EN 15527:2008). It has been tested for ruggedness.

This document is applicable and validated for several types of matrices as indicated in Table 1 (see also Annex A for the results of the validation).

Table 1 — Matrices for which this European Standard is applicable and validated

Matrix	Materials used for validation	
Soil	Sandy soil Mix of soil from the vicinity of Berlin, Germany and PCB-free German reference soil	
Sludge	Mix of municipal waste water treatment plant sludge from North Rhine Westphalia, Germany	
Biowaste	Mix of compost from the vicinity of Berlin, Germany and sludge from North Rhine Westphalia, Germany	
Waste	Contaminated soil, building debris, waste wood, sealant waste, electronic waste, shredder light fraction, cable shredder waste	

WARNING — Persons using this document should be familiar with usual laboratory practice. This document does not purport to address all of the safety problems, if any, associated with its use. It is the responsibility of the user to establish appropriate safety and health practices and to ensure compliance with any national regulatory conditions.

IMPORTANT — It is absolutely essential that tests conducted according to this document be carried out by suitably trained staff.

1 Scope

This document specifies methods for quantitative determination of seven selected polychlorinated biphenyls (PCB28, PCB52, PCB101, PCB118, PCB138, PCB153 and PCB180) in soil, sludge, sediment, treated biowaste, and waste using GC-MS and GC-ECD (see Table 2).

Target analyte CAS-RN a PCB28 2,4,4'-trichlorobiphenyl 7012-37-5 2,2',5,5'-tetrachlorobiphenyl PCB52 35693-99-3 **PCB101** 2,2',4,5,5'-pentachlorobiphenyl 37680-73-2 **PCB118** 2,3',4,4',5-pentachlorobiphenyl 31508-00-6 **PCB138** 2,2',3,4,4',5'-hexachlorobiphenyl 35065-28-2 **PCB153** 2,2',4,4',5,5'-hexachlorobiphenyl 35065-27-1 **PCB180** 2,2',3,4,4',5,5'-heptachlorobiphenyl 35065-29-3 CAS-RN Chemical Abstracts Service Registry Number.

Table 2 — Target analytes of this European Standard

The limit of detection depends on the determinants, the equipment used, the quality of chemicals used for the extraction of the sample and the clean-up of the extract.

Under the conditions specified in this document, lower limit of application from $1 \mu g/kg$ (expressed as dry matter) for soils, sludge and biowaste to $10 \mu g/kg$ (expressed as dry matter) for solid waste can be achieved. For some specific samples the limit of $10 \mu g/kg$ cannot be reached.

Sludge, waste and treated biowaste may differ in properties, as well as in the expected contamination levels of PCB and presence of interfering substances. These differences make it impossible to describe one general procedure. This document contains decision tables based on the properties of the sample and the extraction and clean-up procedure to be used.

NOTE The analysis of PCB in insulating liquids, petroleum products, used oils and aqueous samples is referred to in EN 61619, EN 12766-1 and EN ISO 6468 respectively.

The method can be applied to the analysis of other PCB congeners not specified in the scope, provided suitability is proven by proper in-house validation experiments.

2 Normative references

The following documents are referred to in the text in such a way that some or all of their content constitutes requirements of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN ISO 5667-15, Water quality — Sampling — Part 15: Guidance on the preservation and handling of sludge and sediment samples (ISO 5667-15)

EN ISO 16720, Soil quality — Pretreatment of samples by freeze-drying for subsequent analysis (ISO 16720)

EN ISO 22892, Soil quality — Guidelines for the identification of target compounds by gas chromatography and mass spectrometry (ISO 22892)

ISO 8466-1, Water quality — Calibration and evaluation of analytical methods and estimation of performance characteristics — Part 1: Statistical evaluation of the linear calibration function



This is a free preview	 Purchase the entire 	e publication at the link below:
------------------------	---	----------------------------------

Product Page

- Dooking for additional Standards? Visit Intertek Inform Infostore
- Dearn about LexConnect, All Jurisdictions, Standards referenced in Australian legislation